



The suction liner lid has the **large handle (1)** for easy take it off from the canister and move the liner bag to waste. It has the **large pour spout with 35 mm. diameter (2)** for easy empty without lid remove & for specimen collector. **The patient port (3)** allows direct patient tube connection via **elbow connector (7)** allow quick hook up and prevent kinked tubing. The **tandem port (4)** for multiple suction bag connection. The **vacuum port (5)** for connect the vacuum source .The **closing plugs (6)** for opening ports need closed.

Double safety design to ensure no liquid is overflow to the suction unit and patient tubing ,effectively prevent cross infection and protect the medical staff safe. All FilTrex suction liner bags are equipped with **PTFE filter shut off valve (1)** prevents overflow and bacterial to suction unit and **Check Valve (2)** prevents overflow to patient tube.

Item No.	Suction Liner Bag : Size & Specs	Packaging
SB-1200-01	1,200 ml with Shut-off Valve	50 pcs / pack
SB-1200-02	1,200 ml with Shut-off Valve & check valve	50 pcs / pack
SB-2200-01	2,200 ml with Shut-off Valve	50 pcs / pack
SB-2200-02	2,200 ml with Shut-off Valve & check valve	50 pcs / pack
SB-3200-01	3,200 ml with Shut-off Valve	50 pcs / pack
SB-3200-02	3,200 ml with Shut-off Valve & check valve	50 pcs / pack
Item No.	Accessories : Size & Specs	Packaging
F-EC-01	Elbow connector	12 pcs / pack
CJ-01/02/03	Canister for bag 1.2 liter, 2.2 liter and 3.2 liter	1 pcs / pack
F-WB-01	Wall bracket ,for wall mounting	12 pcs / pack
F-SST-01	Sterile connecting tube with rib for Patient port ,ID 6 mm. x L1.8 m.	25 pcs / pack
F-SST-02	Sterile connecting tube with rib for Patient port ,ID 6 mm. x L 3 m.	25 pcs / pack
F-STT-01	Sterile connecting tube for Tandem port ,ID 6 mm. x L 0.4 m.	25 pcs / pack
F-MSS-01	Mobile suction stand	1 pcs / pack
F-RP-01	Regulator & Pipe connector	1 pcs / pack

NOTE : We reserve the right to alter the specifications of the products described in this leaflet without prior notice.